



CareSens II Meter



i-SENS, Inc.
465-6 Wolgye-dong, Nowon-gu
Seoul 139-845, Korea
www.i-SENS.com



Medical Technology Promedt
Consulting GmbH
Altenhofstrasse 80 D-66386
St. Ingbert, Germany



© 2009 i-SENS, Inc. All rights reserved. PGA1E3013 REV0

Quick Reference Guide

CareSens™ II

Blood Glucose Monitoring System

Accurate, complete blood glucose monitoring through advanced biosensor technology



This guide offers at-a-glance information about how to use your CareSens II Blood Glucose Monitoring System. It's recommended that you read the entire *CareSens II User Manual* before using the system.

Important Safety Information

When using the CareSens II Blood Glucose Monitoring System, these safety rules *must* be followed to prevent injury or infection.


- Before and after getting blood samples or using the system, wash your hands and the intended sample site with soap and warm water. Rinse and dry thoroughly.
- Do *not* use the system for the diagnosis of diabetes without the guidance of a healthcare professional.
- Use *only* fresh capillary whole blood samples for testing.
- A lancet should be used *only* one time. Do *not* use a used lancet.
- Always use a new sterile lancet and keep the lancing device clean.

Coding Your Meter


Please read carefully before using the meter.

1. Gently push test strip in until the meter beeps.



2. See if code numbers on the display and test strip bottle match. If so, wait until the  symbol appears on the screen to begin the test. If not, follow the next step.



3. Use the C or M Button to scroll code numbers up or down. When the  symbol appears, the meter is ready.




Blood Glucose Testing

Before testing, wash your hands and the intended sample site with soap and warm water. Rinse and dry thoroughly.

1. Gently push test strip in until the meter beeps.



2. When the  symbol appears on the display, apply the blood sample to the narrow top edge of test strip until the meter beeps.



3. Test result appears after five(5) seconds.



Using the Meter Memory

1. Press the M Button to turn on the meter. The current date and time will appear at the bottom of the screen for two(2) seconds, followed by the 14-day average and number of tests performed.



2. Press the C Button. The most recent test result, its date, and time will appear. Each time you press the C Button, you'll see previous test results, dates, and times displayed in chronological order. To see test results displayed more quickly, press and hold the C Button.



3. Press the M Button to turn the meter off. The meter will turn itself off after two(2) minutes if the M Button is not pressed.



Display and Messages

1. Each time you turn on the meter by inserting a test strip or pressing the M Button, all segments of the display will briefly appear.



2. If the test result is above 600 mg/dL(33.3 mmol/L), Hi appears on the display.





3. If the test result is below 20 mg/dL(1.1 mmol/L), Lo will appear on the display.



Display and Messages

4. The test strip has been used. Use a new strip.



5. The blood sample was applied before the  symbol appeared. Repeat the test with a new strip. Wait until the  symbol appears before applying blood sample.



6. The temperature was outside the recommended range. Take system where the temperature is within operating range(10-40°C/50-104°F) and repeat the test after 30 minutes.



7. The blood sample did not fill the confirmation window of the test strip during measurement because of abnormally high viscosity or insufficient volume. Repeat the test with a new test strip and adequate blood sample.



Contact your authorized i-SENS sales representative if error messages persist after testing several times.